

Melting Point Apparatus MPA100 — Automated melting point system



1EA

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer [MASSPECS](#)

Description

Overview

OptiMelt provides a fast and accurate means of automatically determining the melting points and melting ranges of chemical substances. With microprocessor-controlled temperature ramping, a built-in digital camera, and a selling price that is half that of competing models, OptiMelt offers the best value of any commercially available melting point apparatus.

MPA100 Specifications

Operation

Temperature Melting point and melting point range

Temperature

display

Start (Ambient + 10 °C) to 396 °C, 0 °C min. set
temp. point

range

Stop (Start temperature + 4 °C) to 400 °C

temp.

range

Temperature 0.1 °C

Temperature

resolution

Ramp 0.1 °C to 20 °C per minute (0.1 °C increments)

rate

Heat-up 10 minutes (50 °C to 350 °C)

time

Cool- 10 minutes (350 °C to 50 °C)

down

time

Temperat±0.3 °C (up to 100 °C)

ure ±0.5 °C (up to 250 °C)

accuracy ±0.8 °C (up to 400 °C)

Reproduc0.2 °C

ibility

TemperatPt RTD (built-in)

ure

sensor

Oven Closed-loop PID

control

General

Display Back-lit, touchscreen LCD (5.8")

Printer RS-232 serial port. Supports Epson compatible

interface dot-matrix printers.

Compute USB. All instrument functions can be queried

r and controlled through a high-level command

interface set.

Capillaries

Dimens 1.4 mm to 2.0 mm outside dia., 100 mm length

ions

Capacit Up to 3 tubes simultaneously

y

Fill 2 mm to 3 mm

height

Power 90 to 264 VAC, 47 to 63 Hz, 125 W

Operatin 0 °C to 40 °C, non-condensing

g temper

ature

Weight 9 lbs.

Dimensio 7.5" x 10" x 8.5" (WHL)

ns

Warranty One year parts and labor on defects in

materials and workmanship

MPA100 Melting Point Apparatus

OptiMelt provides a fast and accurate means of automatically determining the melting points and melting ranges of chemical substances. With microprocessor-controlled temperature ramping, a built-in digital camera, and a selling price that is half that of competing models, OptiMelt offers the best value of any commercially available melting point apparatus.

Automated Measurements

OptiMelt is specifically designed for unattended operation. It has a built-in digital camera that continuously captures real-time images of the samples, and it uses digital image processing to determine results. The melting points and melting ranges are prominently displayed on the front panel and automatically recorded into memory for later review.

**Set
up
Me
nu**

Opt ion Me nu

Simple Operation

OptiMelt has an intuitive front panel and is very easy to use. Simply select a start temperature, ramp rate, stop temperature, and hit start. Results can be easily seen from across the lab on the large LCD display. Samples can be viewed on the front panel through a removable magnification lens. During a measurement, you can flag relevant events by pressing dedicated front-panel buttons. Up to six individual temperatures can be tagged for each sample.

Interactive help is available for all functions and parameters of the instrument. Text and numerical entry keypads are built into the touchscreen interface so that no external keyboard is required.

Fast Measurements

The small, aluminum oven design, along with microprocessor-controlled temperature ramping, provides fast and repeatable warm-up and cool-down cycling. Programmable ramp rates from 0.1 °C/min to 20 °C/min, in 0.1 °C/min increments, provide measurement flexibility.

The ability to rapidly preheat the oven to a start temperature slightly below an expected melting point minimizes analysis time.

Accurate Results

OptiMelt uses a platinum RTD sensor and makes temperature measurements to 400 °C with 0.1°C resolution. It is easily calibrated in the field against certified reference standards and complies date of the last calibration which is included in all reports.

Data Logging

Opt ion O1 00P Pri nte r

OptiMelt meets modern GLP requirements for melting point determination. Up to 24 different analysis methods can be stored in memory. Eight complete melting point reports can also be saved. Records can be displayed on the front panel, printed, or transferred to a host PC via USB. Printouts include technician identification and space for a signature.

Digital Image Processing

The OptiMelt system is unique among melting point analyzers. It has a built-in, high-resolution digital camera to capture real-time images of the samples. It then uses digital image processing to determine phase transitions in the samples.

Unattended melting points, and melting point ranges determined by OptiMelt, closely match visual results. OptiMelt's digital image processing algorithm provides a dramatic improvement over optical absorption or reflection techniques.

Fig

ure
1 -
Solid
Sample
s

Figure
2 -
Meniscus
Point

Figure
3 -
Clear
Point

MeltView Software

The MeltView Window software displays live, high-resolution images of each sample during analysis.

The three screen-captured images to the right were collected during a typical melt. Figure 1 shows the unmelted solid samples, figure 2 reveals the meniscus point, and figure 3 corresponds to the end of the melt (clear point). Stored images may be recalled at any time, and can be played back frame-by-frame or as a movie.

SRS also offers a paid FDA 21 CFR Part 11 ready version. See separate tab above for details.

Thermometer Input and Timersp

OptiMelt is the only melting point apparatus that can also be used as a general-purpose digital thermometer and timer. In between melting point determinations, OptiMelt may be used with a Pt RTD to measure and log temperatures around the lab. Three built-in timers are also available to keep track of timed events during experiments.

Reviews

There are yet no reviews for this product.