

Collagen Gel Culturing Kit



[Kit Contents] (1) Cellmatrix Type 1-A 1 bottle x 20 mL (2) 10 x conc. Ham's F12 1 bottle x 5 mL (3) 10 x conc. MEM 1 bottle x 5 mL (4) Reconstitution buffer 5 bottles x 4 mL

Rating: Not Rated Yet

Price

Sales price 1960 ?

Discount

[Ask a question about this product](#)

Manufacturer [FUJIFILM Wako](#)

Description

This product is for research use only. Do not administer it to human.

Cellmatrix was developed as a collagen for tissue culture, and makes it possible to culture various collagen.

Collagen Gel Culturing Kit (#638-00781 (1 kit))

for Cultivation on Collagen Gel and for Collagen Gel Embedded Culture Method

[Features]

? recommended to beginners for collagen gel culture

? sterile

[Kit Contents]

(1) Cellmatrix Type 1-A 1 bottle x 20 mL

(2) 10 x conc. Ham's F12 1 bottle x 5 mL

(3) 10 x conc. MEM 1 bottle x 5 mL

(4) Reconstitution buffer 5 bottles x 4 mL

[Related Products]

? **Cellmatrix Type 1-P** (#638-00661 (20 mL); #634-00663 (100 mL))

for Cultivation on Collagen Gel, for Collagen Gel Embedded Culture Method and for cultivation on Collagen Coated Dishes

[Features]

? pepsine-solubilized type 1 collagen solution (3 mg/mL, pH 3.0)

? derived from porcine tendon

? gives transparent gels with consistent high gel strength after reconstitution

? **Cellmatrix Type 1-A** (#631-00651 (20 mL); #637-00653 (100 mL))

for Cultivation on Collagen Gel and for Collagen Gel Embedded Culture Method

[Features]

? acid-soluble type 1 collagen solution (3 mg/mL, pH 3.0)

? derived from porcine tendon

? gives transparent gels with consistent high gel strength after reconstitution

? sterile

Cellmatrix LINE-UP, manufactured by Nitta Gelatin Inc. (Osaka, Japan)
Shipped at 2 ~ 10 degrees C with cooling medium.

Reviews

There are yet no reviews for this product.