

## DL-3-Aminoisobutyric Acid Monohydrate Reference Material

The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.

500 mg

Rating: Not Rated Yet

**Price**

Sales price 5050 ?

Discount

*On Order*

★★★★★

[Ask a question about this product](#)

Manufacturer [FUJIFILM Wako](#)

### Reviews

There are yet no reviews for this product.