

Activator Solution-3 (0.25mol/l 5-Benzylthio-1H-tetrazole, Acetonitrile Solution)

The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.

500 ml

Rating: Not Rated Yet

Price

Sales price 950 ?

Discount

[Ask a question about this product](#)

Manufacturer [FUJIFILM Wako](#)

Reviews

There are yet no reviews for this product.