

ANIT HUMAN BFGFM.A. CL#3H3 IGG WHOLE PEROXIDASE CONJUGATED

The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.

200ug

Rating: Not Rated Yet

Price

Sales price 3480 ?

Discount

On Order

★★★★★

[Ask a question about this product](#)

Manufacturer [FUJIFILM Wako](#)

Reviews

There are yet no reviews for this product.