

Poly-D-lysine Hydrobromide MW 30,000 - 70,000 (?-ray irradiated for sterilization)

The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.The Wako logo consists of the word "Wako" in white, sans-serif font, centered within a solid red square.

5 mg

Rating: Not Rated Yet

Price

Sales price 670 ?

Discount

On Order



[Ask a question about this product](#)

Manufacturer [FUJIFILM Wako](#)

Reviews

There are yet no reviews for this product.